

0.005N Sulfuric Acid Mobile Phase

Ready-to-Use Solutions | Acidified Mobile Phases

1. Specification scope

Product	0.005N Sulfuric Acid Mobile Phase
CAS number	7664-93-9
Grade scope	HPLC Prep / HPLC Gradient / HPLC Isocratic
Document status	Reference specification summary

The values below are representative catalog specifications. Where several grades are listed, acceptance limits may differ by grade. The current controlled grade-specific specification agreed at quotation or qualification is the governing document.

2. Test requirements

Test item	Specification / reference value	Test method
H ₂ SO ₄ Concentration	0.005N (±2%)	Titration
Water Quality	TOC < 50 ppb, Resistivity > 18 MΩ·cm	TOC/Conductivity
pH	2.2 ± 0.1	pH meter
Filtration	0.2 μm membrane filtered	Visual
Packaging Atmosphere	Nitrogen-filled	-
Compatible Columns	Aminex HPX-87H, Bio-Rad, Shodex SH-1011	-
Detects	Acetic acid, lactic acid, citric acid, succinic acid, malic acid	Approved method
Column Type	C18 reversed-phase	Approved method
Grade/Performance	Ready to use , Sterile filtered , pHStable , Low impurities	Specification
Product Scope	H ₂ SO ₄ /Ultrapure water , ~2.0	Catalog
Available Packaging	5L, 7L, 10L, 20LAluminum composite bag(Nitrogen filled)	Catalog
Documentation	CoA/SDS/TDS available on request	Document control

3. Packaging and release documentation

Package size	Container / format	Supply status
5L	FLEX BAG	Available
7L	FLEX BAG	Available
10L	FLEX BAG	Available
20L	FLEX BAG	Available

Each released batch is supported by a batch-specific Certificate of Analysis. Packaging sizes are consolidated in this document; actual availability, closure system and minimum order quantity are confirmed with the quotation.

4. Document control notice

This downloadable file is a reference summary generated from the LANCHROM product catalog. Request the current approved specification, SDS and sample COA before method validation, supplier qualification, regulatory submission or purchase release.